

News Issue:

Modular Aseptic Clean Room Facility for Tameside General

Mach-Aire Ltd have designed and installed on site at Tameside General Hospital a new Pharmacy Department Aseptic Clean Room Facility. The new building has been constructed adjacent to the Pharmacy Department and provides state of the art fully GMP compliant facilities for the dispensing of sterile medical products.

The facility was constructed at Mach-Aire's factory and delivered in three modular sections onto pre-formed foundations. Included were two Grade 'B' isolator rooms each fitted with a Mach-Aire 4-glove positive pressure isolator.

The facility comprises, Outer Support, Change Rooms, Inner Support, CIVAS and TPN Isolator Rooms, all terminal HEPA filtered and air conditioned.

Innovations such as timed transfer hatches, alcohol extract stations and a sophisticated CCTV system combined with a clean room intercom system allow for good staff working conditions with a high level of supervision and security for the staff.



Fume Cupboards & Extraction Systems, Pharmacy Aseptic Clean Rooms, Isolators and Powder Containment