

## Technical Specifications – Laminar Flow Cabinets

### Horizontal Laminar Flow

MODEL	Overall DIM (LxDxH mm)	WORKSPACE DIM (LxDxH mm)	AIR VOL EXHAUST	SHIPPING WEIGHT
HFC-900	1000 X 900 X 1180	995 X 595 X 770	0.243m <sup>3</sup> /sec	115 kg
HFC-1200	1300 X 900 X 1180	1195 X 595 X 770	0.324m <sup>3</sup> /sec	150 kg
HFC-1500	1600 X 900 X 1180	1495 X 595 X 770	0.405m <sup>3</sup> /sec	205 kg
HFC-1800	1900 X 900 X 1180	1795 X 595 X 770	0.486m <sup>3</sup> /sec	255 kg

