



MACH-AIRE

WEIGH STATION



The weigh station is a neat solution for the containment of fine powders that may be released whilst dispensing/weighing during normal laboratory procedures.

FEATURES

- Excellent containment
- No effect on weigh scales readings up to 5 decimal places
- Low flow
- Clean surfaces
- Real time performance monitoring.


APPLICATIONS

- Pharmaceutical research
- Chemistry
- Healthcare
- Toxic powder weighing.

DIMENSIONS

MODEL	WIDTH	EXTERNAL FLOW
600WS	600mm	0.056m ³ /s
900WS	900mm	0.840m ³ /s
1200WS	1200mm	0.112m ³ /s

Available in filtered and with in built recirculation fan.



INNOVAIRTION®

Mach-Aire Limited
 Bridge Street
 Horwich
 Bolton BL6 7BT

Tel: +44 (0) 1204 668905
 Fax: +44 (0) 1204 668906
 e-mail: sales@machaire.co.uk
www.machaire.co.uk

Registered in Scotland 57500 VAT Registered 592 9715 90